

## **Technical Data Sheet**

Farapol Jam Chemical Indus. Co.

## FARAPOL O 137

## **NPG Resin**

Product DescriptionFarapol O 137 is an Unsaturated Polyester Resin based on Neopentyl Glycol and Orthophthalic<br/>anhydride, dissolved in and cross-linked with Styrene Monomer. It is medium reactive, non-accelerated<br/>and non-thixotropic resin. This resin can be supplied pre-accelerated, with a special type of colorless<br/>Cobalt Octoate

Applications and<br/>UseThis resin is designed for manufacturing of artificial marble, architectural panels, filled resin objects,<br/>decorative castings, vanities, bathtubs and wall panels.

**Certificates and** Approvals Farapol O 137 is synthesized from raw materials listed in FDA regulation Title 21 CFR 177.2420. Farapol Jam Chemical Industrial Company carries out this resin's production, quality control, and distribution in compliance with ISO 9001, 14001, 45001, 10002, 10004, 10015, and 17065 standards.

| Typical Liquid   | Property @ 25 °C                  | Unit     | Specification | Method      |
|------------------|-----------------------------------|----------|---------------|-------------|
| Resin Properties | Viscosity Brookfield <sup>1</sup> | cps      | 500 - 550     | ISO 2555    |
|                  | Acid Value                        | mgKOH/g  | Max 28        | ASTM D 1639 |
|                  | Solid Content                     | %        | 62 - 64       | ISO 3251    |
|                  | Color                             | Gardner  | Max 1         | ASTM D 1544 |
|                  | Specific Gravity                  | relative | 1.12-1.14     | ISO 2811    |
|                  | Gel Time <sup>1</sup>             | minute   | 8 - 10        | ASTM D 2471 |
|                  | Exothermic Peak Temperature       | °C       | 150-180       | ASTM D 2471 |

1) Gel Time and Viscosity can be adjusted as per customer requirements.

| Dehavior of   |  |
|---|--|
| Benavior of<br>Resin <sup>2</sup> Gel Time (minute)20-258-105-7 |  |

2) Mix ratio for measuring Gel Time: (Cobalt Octoate Farapol C 901 1%- 1.0 phr, Akperox A60 1.0 phr).



| Typical Casted   | Property   | Unit           | Specification       | Method                              |  |  |  |
|--|--|----------------|---------------------|-------------------------------------|--|--|--|
| Resin Properties <sup>3</sup>                            | Tensile Strength   | MPa            | Min 70              | ISO 3268, ASTM D638,<br>ISO 527-2&4 |  |  |  |
|  | Elongation at Break  | %              | Min 3.0             | ISO 3268, ASTM D638,<br>ISO 527-2&4 |  |  |  |
|  | Tensile Modulus  | GPa            | Min 3.0             | ISO 3268, ASTM D638,<br>ISO 527-2&4 |  |  |  |
|  | Flexural Strength  | MPa            | Min 120             | ISO 178/ASTM D 790                  |  |  |  |
|  | Flexural Modulus   | GPa            | Min 3.5             | ISO 178/ASTM D 790                  |  |  |  |
|  | Heat Deflection Temperature (HDT)  | <sup>0</sup> C | Min 70              | ISO 75-2 Test Method A              |  |  |  |
|  | Barcol Hardness  | Barcol         | Min 40              | ASTM D 2583                         |  |  |  |
|  | Water Absorption   | %              | pprox 0.20          | ISO 62- Test Method 3               |  |  |  |
|  | Linear Shrinkage <sup>4</sup>  | %              | $\approx 1.6$       | Internal method                     |  |  |  |
|  | <ul> <li>3) Materials used for curing are: (Cobalt Octoate Farapol C 901 1% - 1.0 phr, Akperox A60 1.0 phr). Curing Time is 2 Room Temperature and 3 hrs at 80 °C.</li> <li>4) This test has been done on the specimen with linear dimensions (1 cm × 1 cm × 100 cm).</li> </ul>   |                |                     |                                     |  |  |  |
| Handling, Storage<br>and Stability<br>Healthy and Safety | <ul> <li>FARAPOL O 137 is a product that is sensitive to temperature, light, and oxidation. Hence, it should be stored indoors in a dry place at a temperature between 5 and 25°C. Keep always in the original, unopened, and undamaged containers. Avoid keeping material exposed to sunlight. On storage under the above-mentioned conditions, the shelf life for FARAPOL O 137 is 6 months.</li> <li>Avoid storing the resin along with Metallic Driers and Peroxides in the same area. Safety Datasheets of the product are available on demand. The user is fully responsible for reviewing the material's Safety Data Sheet (SDS) and understanding proper handling procedures prior to using the product.</li> </ul>  |                |                     |                                     |  |  |  |
| Packaging  | FARAPOL O 137 is supplied in 200   | Kg steel bar   | rels, IBC tanks and | bulk road tankers.                  |  |  |  |
| 8 8  |  | U              | ,<br>,              |                                     |  |  |  |
| Notice   | The information contained herein is provided in good faith and is with the best of our accurate knowledge, but we assume no liability for its accuracy or completeness. Therefore, the buyer is advised to determine the suitability of this product for the intended use. We retain the right to make any changes according to technological progress or further developments. Please note that variations in testing conditions across different laboratories may result in discrepancies, and a tolerance of up to 5% in test results should be expected. Farapol Jam Company reserves the right to modify the information in this document at its discretion. The latest version available on the Farapol website is considered valid, and any previous versions are void. |                |                     |                                     |  |  |  |
| Document<br>Registration                                 | Pub. No: POL- F-76-33 Revisio  | on No.: 3      | Rev. Date: 01/1     | 8/2025                              |  |  |  |
| Contact Information                                      | Tel: +98 21 26231019 Fax   | : +98 21 26    | 231014              |                                     |  |  |  |

Site: www.farapol.com

Email:info@farapol.com

